



## POCKET BALANCE (0.01g) 100012

A pocket sized balance with 0.01g resolution. The entire balance is contained within a protective hard folding case for travel. During operation the cover folds to reveal instructions conveniently printed on the inside, a stainless steel weighing platform, and a weighing dish. Features five units of measure, automatic calibration, automatic shut off, and a tare function. The display indicates overload and low battery and freezes automatically when stable for 2 seconds. If weight is added or subtracted Pocket Balance resumes measuring.

Dim: 5½" x 3 x ¾" (138 x 75 x 20mm), platform 2¾" x 2⅜" (70 x 60mm)

Weight: 7oz (196g) with 4 AAA batteries (included).

## OPTIONAL ACCESSORIES

100013 100 gram calibration weight

Units	Range	Resolution
Grams	100	0.01
Ounces	3.5	0.001
Pennyweight	64	0.01
Troy Ounces	3	0.001
Carat	500	0.05