



MINIMIXER 460001

Stir in test tubes, micro centrifuge tubes, and other containers too narrow for conventional mixers. The stirring rod will also flare out for use in larger spaces. Clear vinyl carrying case, 4½ x ¼" stirring rod. Batteries not included. Weight: 3oz (85g). Dimensions (including rod): 8" x 1¼" x 1" (203 x 32 x 25mm).

ADDITIONAL STIRRING RODS 460002

Disposable nylon stirring rods for use with Minimixer.

SILICONE SAFETY PIPET FILLERS

Enables filling and dispensing of pipets with one hand, safely through avoidance of hand or mouth contact with liquid. Our Silicone Pipet Fillers offer a softer feel with easier and faster operation of valves. Silicone is also more chemically resistant and more resilient than either PVC or natural rubber. It will not crack or become brittle and will last longer and stay flexible. Adjusts to fit any size pipet. Chemically resistant glass valves. Available in two colors for identification purposes (470113-Red, 470114-Black). Glass tube not included. Weight: 1oz (37g). Dimensions: 4½" x 2" x 2" (114 x 51 x 51mm).

470113 Silicone Pipet Filler, Red

470114 Silicone Pipet Filler, Black

UNIVERSAL BLACK RUBBER SAFETY PIPET FILLER 470221