



### WIRELESS THERMO/HYGROMETER SET, 800029

Read the temperature and RH inside a humidity chamber or desiccator without opening the door. Mount the sensor in green houses, warehouses, air ducts and other hard to reach places and continuously monitor their ambient temperature and RH from a remote location. The sensor broadcasts temperature & RH to the monitor from up to 80' (25M) away (in an open area) via a wireless RF signal. The monitor and sensor each simultaneously display the time in large digits, RH and temperature in small digits. The monitor shows both local and remote readings. Reads temperature in °C and °F, has an alarm clock function and indicates low battery in both sensor and monitor. Both the monitor and sensor can be wall mounted and the monitor also comes with a stand for bench top use. Batteries are included. Weight, monitor: 15oz (427g), sensor: 3oz (91g). Dim, monitor: 10" x 7½" x 1" (255 x 190 x 30mm), sensor: 4¼" x 2½" x 1" (107 x 65 x 30mm).

	Monitor	Sensor
<b>RH Range:</b>		20~90%
<b>Resolution:</b>		1%RH
<b>Accuracy:</b>		5%RH
<b>Temperature Range:</b>	0~50°C & 32~122°F	-20~60°C & -4~140°F
<b>Resolution:</b>	0.1°	0.1°
<b>Accuracy:</b>	1°C, 2°F	1°C, 2°F
<b>Batteries:</b>	3 X AA	2 X AAA



With NIST traceable certificate of compliance 800029C

