



MINI TIME CONTROLLER 810030

Especially useful with recorders, electrophoresis equipment, pumps, and for control and safety on equipment with heating elements such as hot plates, water baths, ovens, and sterilizers.

Plug any piece of laboratory equipment into Time Controller and you can now program it for up to 6 on and 6 off settings per day, a total of 84 settings a week. The programs are user friendly and easy to set with the aid of convenient buttons on the front panel. They can be set to run for as little as one minute in duration. Also offers manual override. An LED light indicates that power has been switched on. The LCD displays real time as well as programmed times. Precise and stable timing control is guaranteed by a crystal oscillator. In case of a power failure the back-up battery (1.5V UM-5 size) ensures that your programs will not be lost (included).

Time Controller will operate any device with a maximum load of 15 Amp, 1800 Watt, and any motor up to 1 HP. Comes with easy to follow instructions. Dimensions: 4" x 2" x 3¼" (102 x 81 x 83mm). Weight: 5oz (142g).



With NIST traceable certificate of compliance 810030C

