



BENCH-TOP PH / MV METER 860031

Rugged housing with large LCD displaying the parameter being measured together with time, date and temperature (in °C or °F).

Features:

- Automatic or manual temperature compensation
- Automatic or manual ranging
- Automatic buffer recognition
- Digital and analog outputs.
- Min/Max and 99 data points
- Data Hold
- 5 calibration points
- Electrode status indication

Operates on AC voltage (no batteries). Comes ready to use with ATC pH electrode, electrode holder arm, AC adaptor, computer cable & software. Also accepts any pH or ORP probe with a standard BNC connector. Weight: 18oz (533g).

Dimensions: 8½" x 6½" x 2¼" (217 x 168 x 58mm).

Unit of Measure	Range	Resolution	Accuracy
pH	0~14	0.01	±0.02
MV / ORP	-1999~+1999	0.1mv -199.9~+199.9 otherwise 1mv	±0.2mv -199.9~+199.9 otherwise ±2mv
Temperature °F	23~176	0.1	±0.9 °F
Temperature °C	-5~80	0.1	±0.5 °C

OPTIONAL ACCESSORIES

- 840016 pH Probe (non ATC)
- 840049 Spear Tip pH Probe
- 840094 USB Cable
- 850059P Replacement ATC pH Probe
- 850088 ORP Probe
- 860008 pH 4, 3 bottles 40 ml each
- 860009 pH 7, 3 bottles 40 ml each
- 860010 pH 10, 3 bottles 40 ml each
- 860011 Deionized Water, 3 bottles, 40 ml each

