



DIAL HYGROMETER/THERMOMETER 736920

For classroom, lab or industrial use. Large, clear, 20 ~ 100% RH scale on 4" dial. Easily read from a distance. Uses a precise synthetic hair hygrometer with an accuracy of $\pm 5\%$. Plastic face meets food industry requirements. Inner dial shows temperature using a bimetal thermometer, with an accuracy of $\pm 2^\circ\text{F}$. Wall mounts. Weight: 3 $\frac{3}{4}$ oz (106g).



With NIST traceable certificate of compliance 736920C

DIAL BAROMETER 736930

Large, clear dual-scale in hPa (hecto Pascal) and mm (millimeter mercury) on a 4" dial which is easy to read from a distance. For classroom, lab, or industrial use. Accuracy of ± 3 to 4 hPa. Wall mounts. Weight: 2 $\frac{2}{3}$ oz (74g).

